

LibraryDoser®

Fully Automated Powder Dispensing

- Walkaway, Fast, Gravimetric Dosing!
- User specifies target masses
- 20 Changeable Sources
- Target Labware = both Minitubes & Vials
- 4 place Analytical Balance standard
- Sub-milligram Precision
- Small Form Factor, 18"x18"x18" configuration dependent.
- Fits in Hood, Glovebox or Enclosure
- Imports/Exports using CSV files
- Comprehensive Run Data for Every Dispense
- 21CFR11 Compliance Capable
- Emails User Upon Completion
- Standard Labware & Inexpensive Consumables



siriusautomation.com/LibraryDoser – for more information

Sirius Automation Group LLC, 318 W. Half Day Road, Buffalo Grove, IL 60089, USA

Tel 847.607.9378

Fax 847.607.9369

www.siriusautomation.com